



Compact Weekday Pill Dispenser Usage and Maintenance Instructions

Ver. 1 Dec. 2011
2223

VM930D **Pink Compact Weekday Pill Dispenser**

VM930E **Blue Compact Weekday Pill Dispenser**

The handy, compact Weekday Pill Dispenser provides a simple reminder to take your medication, with one clearly-labelled compartment per day. It is ideal for use at home, during travel or holidays and will fit neatly into a pocket or bag.

1. Introduction

Thank you for deciding to purchase the Compact Weekday Pill Dispenser from Aidapt. The Compact Weekday Pill Dispenser is manufactured from the finest materials available. When used correctly it is designed to give many years of reliable, trouble-free service.

2. Intended Use

The handy, Compact Weekday Pill Dispenser provides a simple reminder to take your medication, with one clearly-labelled compartment per day.

3. Using your Pill Box Keyring

The round Compact Weekday Pill Dispenser has a simple-to-use operation:

1. Click on the button to rotate the compartments around and select the correct day
2. Flip down the lid to insert or tip out the tablets
3. Click the lid down into place to close the compartment

4. Cleaning

Clean your Compact Weekday Pill Dispenser using a non-abrasive cleaner or mild detergent with a soft cloth only. Only use a mild detergent or disinfectant; never use an abrasive cleaner. Always ensure that you wipe dry all equipment.

5. Reissue

If you reissue or are about to reissue this product, please thoroughly check all components for their safety.

If in any doubt, please do not issue or use, but immediately contact your supplier for service support.

6. Important Information

The information given in this instruction booklet must not be taken as forming part of or establishing any contractual or other commitment by Aidapt Bathrooms Limited, Aidapt (Wales) Ltd or its agents or its subsidiaries and no warranty or representation concerning the information is given.

Please exercise common sense and do not take any unnecessary risks when using this product; as the user you must accept liability for safety when using the product.

Please do not hesitate to contact the person who issued this product to you or the manufacturer (detailed below) if you have any queries concerning the assembly/use of your product.

7. Care & Maintenance

Please make a safety check of the product at regular intervals or if you have any concerns.